

*POLYSOMNOGRAPHY(PSG) SLEEP STUDY
SAMPLE REPORT*

DIAGNOSTIC SLEEP STUDY REPORT SUMMARY

Patient: DEMO REPORT

AGE:67Yrs/M

Weight:89Kg

Height: 156Cm

DATE:1955-05-11

BMI:34217.6 kg/m2

Indication:


PROCEDURE:

Single overnight unattended diagnostic sleep study (level II) done using portable data acquisition devices.


REPORT:

There were no major technical issues with the study.

The AHI was 67.3 /h with the lowest oxygen saturation recorded was 85.9 % from the baseline of 95%.



Sleepparameter	Noraml value
Stage N1(%TST)	3-8%
Stage N2(%TST)	45-55%
Stage N3(%TST)	15-20%
StageR(%TST)	20-25%



Polysomnography Report

A. Patient Information

PatientID:	20250312-14	Name:	DEMO REPORT	Gender:	Male
Height(cm):	156	Weight(kg):	89	BMI:36.6	34217.6
Neck Circ.(cm):	-	Birth Date:	1955-05-11	Age:	67
Tel.		Start Time:	20:21:00	End Time:	06:56:30
Duration(min):	635.5	ValidDuration(min):	632		


B. Sleep Over view

Start Sleep Time:	20:21:00	Wake Duration(min):	49.5
Total Sleep Time(min):	582.5	Sleep latency(min):	4
REM latency(min):	7	Effect Sleep Percent (%) :	92.2


Sleep Stage Events Form

Stage	Duration(min)	Percent(%)	OSA(min)	CSA(min)	MSA(min)	Hypopnea(min)	OD(min)
WAKE	49.5	7.8	0.4	0.1	-	7.7	30.6
NREM I	336	53.2	53	16.1	-	48.2	240.6
NREM II	59	9.3	7.9	0.9	-	12	35.4
NREM III	24.5	3.9	6.7	2.1	-	1.9	17.3
REM	163	25.8	25.2	16.4	-	19.7	110

C. Sleep Trend Chart



A Complete Healthcare Solution



D. Events Analysis

Arousals Event

Arousals Duration(min):	29.4	Arousals Count:	308
Avg.ArousalsDuration(sec):	5.7	Arousals Index:	30.8
Wake Count:	25		

Breath Event

BreathEvent Duration(min):	210.2	Breath Event Count:	653
Avg. Breath Event Duration(sec):	19.3	AHI:	67.3

1.Apnea Event

Apnea Duration(min):	128.4	Apnea Count:	371
Avg. Apnea Duration(sec):	20.8	Apnea Index(AI):	38.2

Type	Stage	Duration(min)	Count	Avg. (sec)	Longest(sec) / Time	Index
OSA	NREM	67.6	218	18.6	43.9/01:15:47	31.2
	REM	25.2	65	23.2	114.5/01:31:30	23.9
CSA	NREM	19.1	51	22.5	49.9/01:22:46	7.3
	REM	16.4	37	26.7	62.6/04:28:07	13.6
MSA	NREM	-	-	-	-/-	-
	REM	-	-	-	-/-	-
Total	NREM	86.7	269	19.3	49.9/01:22:46	38.5
	REM	41.6	102	24.5	114.5/01:31:30	37.5

2.Hypopnea Event

Hypopnea Duration(min):	81.9	Hypopnea Count:	282
Avg. Hypopnea Duration(sec):	17.4	Hypopnea Index(HI):	29

Event	Stage	Duration(mi)	Count	Avg.	Longest	Happen Time	Index
Hypopnea	NREM	62.2	225	16.6	49.2	21:38:51	32.2
	REM	19.7	57	20.7	84.2	22:58:21	21

SpO2 Event

Max.SpO2:	99	Min.SpO2:	58
Avg. SpO2:	90	Oxygen Desaturation Duration(min):	433.8
<95% Duration(min):	538.1	<95%Percent(%):	85.4
<90% Duration(min):	196.5	<90%Percent(%):	31.2
<85% Duration(min):	81.4	<85%Percent(%):	12.9
<80% Duration(min):	36.4	<80%Percent(%):	5.8
OxygenDesaturationCount:	736	Oxygen Desaturation Index:	73.6

Stage	SpO2 Range	OD	OD Count	Avg.	Longest(sec) / Time
Wake	78-99	30.6	53	34.6	104 / 01:46:08
NREM	58-99	293.2	514	34.2	117 / 22:20:20
REM	61-99	110	169	39.1	115 / 23:01:28

Body Position Information

Position	Duration(min)	Percent(%)	Hypopnea	OSA	CSA	MSA	OD	Total
Supine	213.5	33.8	105	68	18	-	229	420
Prone	1	0.2	-	-	-	-	1	1
Left Side	0.1	0	-	-	-	-	-	-
Right Side	406.2	64.3	177	215	70	-	499	961
Stand	11.1	1.8	-	-	-	-	7	7

Heart BeatInformation

Max.Heart Beat:	115	Min. Heart Beat:	41
Avg. Heart Beat in REM+NREM:	74		

Leg Movement Event

Leg Move Count:	155	Leg Move Index:	15.5
PLMS Count:	6	PLMS Index:	0.6
PLMS+ArousalsCount:	-	PLMS+Arousals Index:	-

Snore Event

SnoreDuration(min):	21.7	Snoring Count:	1081
Avg. SnoreDuration(sec):	1.2	Snoring Index:	111.3

F. Summary Information

Avg. SpO2:	90	Min. Spo2:	58	ODI:	73.6
OSA Count:	283	CSA Count:	88	MSA Count:	-
Apnea Index:	38.2	Hypopnea Count:	282	Hypopnea Index:	29
AHI:	67.3	RERA Count:	42	RDI:	71.6
Snoring Index:	111.3	Avg. Pulse:	73.2	BodyChangeCount:	734

Monitor Summary:

Physician's Signature:

Date:

ApneaLink - Report of 28-02-2025 18:34

Treating physician

Referral to

Patient data


First name:	DEMO REPORT	Patient ID:	27022025
Last Name:	DEMO REPORT	DOB:	28-11-1984
Street:	Godasar	Height:	180 cm
City, ST, Zip:	Ahmedabad 380050	Weight:	80 kg
Phone:		BMI:	24.7 kg/m ²

Recording

Date: 27-02-2025
 Start: 22:54 .
 End: 06:10 .
 Duration: 7 h 17 min

Evaluation

Start: 23:04 .
 End: 06:08 .
 Duration: 4 h 34 min



* See Clinical Guide for abbreviations and ResMed standard parameters

Analysis (Flow evaluation period: 4 h 34 min / SpO₂ evaluation period: 7 h 7 min)


Indices		Normal	Result	
AHI*:	23.2	< 5 / h	Average breaths per minute [bpm]:	10.87
RI*:	26.2	< 5	Breaths:	2977
Apnea index:	10.5	< 5 / h	Apneas:	48
UAI:	2.2		Unclassified apneas:	10 (21%)
OAI:	6.6		Obstructive apneas:	30 (62%)
CAI:	1.5		Central apneas:	7 (15%)
MAI:	0.2		Mixed apneas:	1 (2%)
Hypopnea index:	12.7	< 5 / h	Hypopneas:	58
% Flow lim. Br. without Sn (FL):	37	< Approx. 60	Flow lim. Br. without Sn (FL):	1089
% Flow lim. Br. with Sn (FS):	0	< Approx. 40	Flow lim. Br. with Sn (FS):	5
			Snoring events:	983
ODI Oxygen Desaturation Index*:	18.4	< 5 / h	No. of desaturations:	131
Average saturation:	93	94% - 98%	Saturation <= 90% :	24 min (6%)
Lowest desaturation:	81	-	Saturation <= 85% :	1 min (0%)
Lowest saturation:	81	90% - 98%	Saturation <= 80% :	0 min (0%)
Baseline Saturation:	96	%	Saturation <= 89% :	13 min (3%)
			Saturation <= 88% :	7 min (2%)
Minimum pulse:	52	> 40 bpm		
Maximum pulse:	111	< 90 bpm		
Average pulse:	70	bpm		
Proportion of probable CS epochs:	0	0%		



Analysis status: Analyzed automatically


Analysis parameters used (Default)

Apnea [20%; 10s; 80s; 1.0s; 20%; 60%; 8%]; Hypopnea [70%; 10s; 100s; 1.0s]; Snoring [6.0%; 0.3s; 3.5s; 0.5s]; Desaturation [4.0%]; CSR [0.50]

Comments







RCH PLUS HEALTHCARE PVT. LTD. /+91 9723891128 ECG HOLTER Report

Name: DEMO REPORT		Sex: Male	Age: 66Years
Department: AT HOME	Bed No.: N/A	Admission No.: N/A	Case No.: 20230414
ID: 0000012	Start: 2023-04-14 09:14:08	Stop: 2023-04-16 09:03:26	Total time: 47:49:18

HR: Av. HR: 88 Min HR: 66 Time: 2023-04-15 03:07:16 Max HR: 121 Time: 2023-04-15 03:38:43 Total Beats: 251552 Abnormal Beats: 38011 Abnormal Permillage: 151			Pause time: Pause Beats that R-R > 2000 ms: 0		
--	--	--	--	--	--

S:			V:		
Total S:	477	Single: 468 SVE	Total V:	37534	Single: 23504 VE
Total S Pair:	2	Total: 4 SVE	Total V Pair:	25	Total: 50 VE
Bigeminy:	0	Total: 0 SVE	Bigeminy:	3091	Total: 10441 VE
Trigeminy:	0	Total: 0 SVE	Trigeminy:	924	Total: 3493 VE
Total S Run:	1	Total: 5 SVE	Total V Run:	15	Total: 46 VE
S Permillage:	1		V Permillage:	149	
Max S in a Min:	7		Max V in a Min:	47	

HRV(Time Domain):		ST:	
SDNN (ms) :	105.7		
SDANN (ms) :	57.2		
rMSSD (ms) :	97.1		
PNN50 (%) :	22.6		
CV :	0.12		
HRV(Frequency Domain):(ms*ms/Hz)			
Power:	4592.2		
ULF :	736.8		
VLf :	2104.3		
LF :	845.9		
HF :	905.2		

Conclusion:

Physician Signature:

Arrhythmia Table

Time	HR (bpm)			V (n)					S (n)					L n	BC seg.
	Av.	Min.	Max.	V	Pair	Run	Big	Trig	S	Pair	Run	Big	Trig		
09:14 10:14	93	81	110	506	0	0	52	0	2	0	0	0	0	0	0
10:14 11:14	90	76	109	762	3	0	88	0	7	0	0	0	0	0	0
11:14 12:14	85	76	96	740	1	0	84	0	2	0	0	0	0	0	0
12:14 13:14	88	77	101	428	1	0	26	0	2	0	0	0	0	0	0
13:14 14:14	85	77	95	1113	0	0	106	21	18	0	0	0	0	0	0
14:14 15:14	88	83	95	691	0	0	9	34	34	0	0	0	0	0	0
15:14 16:14	95	85	108	1086	4	0	111	20	25	0	0	0	0	0	0
16:14 17:14	96	83	112	1032	2	2	113	0	5	0	0	0	0	0	0
17:14 18:14	87	80	96	644	0	0	61	0	1	0	0	0	0	0	0
18:14 19:14	86	81	93	250	0	0	15	0	1	0	0	0	0	0	0
19:14 20:14	90	80	101	234	0	0	9	0	2	0	0	0	0	0	0
20:14 21:14	84	79	90	474	1	1	9	0	1	0	0	0	0	0	0
21:14 22:14	90	80	105	543	1	0	43	0	5	0	0	0	0	0	0
22:14 23:14	88	85	91	1189	0	0	144	37	6	0	0	0	0	0	0
23:14 00:14	88	82	96	1194	1	0	86	63	8	0	0	0	0	0	0
00:14 01:14	87	81	94	1168	0	0	76	56	27	0	1	0	0	0	0
01:14 02:14	88	77	101	1250	1	0	18	86	16	0	0	0	0	0	0
02:14 03:14	80	71	90	958	0	0	40	75	28	0	0	0	0	0	0
03:14 04:14	88	71	117	969	1	0	94	29	21	0	0	0	0	0	0
04:14 05:14	80	71	90	923	0	0	47	40	34	0	0	0	0	0	0
05:14 06:14	83	73	95	769	0	0	36	30	28	1	0	0	0	0	0
06:14 07:14	88	71	114	666	0	3	62	5	26	1	0	0	0	0	0
07:14 08:14	94	83	108	793	0	0	114	3	3	0	0	0	0	0	0
08:14 09:14	91	77	112	250	3	0	13	0	9	0	0	0	0	0	0
Current day: 1 / Total 2 days															

Arrhythmia Table

Time	HR (bpm)			V (n)					S (n)					L n	BC seg.
	Av.	Min.	Max.	V	Pair	Run	Big	Trig	S	Pair	Run	Big	Trig		
09:14 10:14	94	81	111	290	0	0	24	0	2	0	0	0	0	0	0
10:14 11:14	89	80	101	102	0	2	1	0	2	0	0	0	0	0	0
11:14 12:14	83	75	93	112	0	0	0	0	0	0	0	0	0	0	0
12:14 13:14	86	74	102	293	0	1	2	0	0	0	0	0	0	0	0
13:14 14:14	88	76	103	213	0	1	11	0	5	0	0	0	0	0	0
14:14 15:14	86	76	97	1172	0	0	109	14	3	0	0	0	0	0	0
15:14 16:14	104	90	122	1105	1	0	92	50	15	0	0	0	0	0	0
16:14 17:14	96	86	107	1073	1	0	151	0	3	0	0	0	0	0	0
17:14 18:14	95	83	110	980	0	0	78	0	1	0	0	0	0	0	0
18:14 19:14	93	82	108	625	0	0	66	0	1	0	0	0	0	0	0
19:14 20:14	92	80	107	291	0	0	17	0	0	0	0	0	0	0	0
20:14 21:14	93	81	111	274	0	1	6	0	1	0	0	0	0	0	0
21:14 22:14	92	81	105	575	1	0	80	0	3	0	0	0	0	0	0
22:14 23:14	89	82	98	1116	1	0	153	22	8	0	0	0	0	0	0
23:14 00:14	93	85	102	1420	0	0	207	1	6	0	0	0	0	0	0
00:14 01:14	84	75	96	884	0	0	94	1	13	0	0	0	0	0	0
01:14 02:14	84	73	98	823	0	0	76	4	21	0	0	0	0	0	0
02:14 03:14	88	83	94	1234	0	0	63	75	8	0	0	0	0	0	0
03:14 04:14	90	80	103	1310	0	0	35	86	15	0	0	0	0	0	0
04:14 05:14	88	75	106	1253	0	0	14	75	17	0	0	0	0	0	0
05:14 06:14	99	84	120	1120	0	1	19	82	31	0	0	0	0	0	0
06:14 07:14	99	84	118	1099	0	3	118	15	4	0	0	0	0	0	0
07:14 08:14	86	73	104	927	1	0	122	0	7	0	0	0	0	0	0
08:14 09:03	88	79	100	611	1	0	97	0	0	0	0	0	0	0	0
All	90	71	122	37534	25	15	3091	924	477	2	1	0	0	0	0
Current day: 2 / Total 2 days															

ST Table

Time	I_ST mV	II_ST mV	III_ST mV	aVR_ST mV	aVL_ST mV	aVF_ST mV	V1_ST mV	V2_ST mV	V3_ST mV	V4_ST mV	V5_ST mV	V6_ST mV
09:14 10:14	0.00	0.00	-0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.00	0.00
10:14 11:14	0.01	0.00	-0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.00	0.00
11:14 12:14	0.01	0.02	0.01	-0.01	-0.00	0.01	0.01	0.03	0.03	0.02	0.01	0.01
12:14 13:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.01
13:14 14:14	0.00	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
14:14 15:14	0.00	-0.00	-0.01	0.00	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
15:14 16:14	0.00	-0.02	-0.02	0.01	0.01	-0.02	0.01	0.02	0.02	0.00	-0.00	-0.01
16:14 17:14	0.00	-0.00	-0.01	0.00	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
17:14 18:14	0.01	0.02	0.01	-0.01	0.00	0.01	0.01	0.03	0.03	0.02	0.01	0.01
18:14 19:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.03	0.02	0.02	0.01	0.01
19:14 20:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.01
20:14 21:14	0.00	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
21:14 22:14	0.01	0.01	0.00	-0.00	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
22:14 23:14	0.00	0.02	0.01	-0.01	-0.00	0.01	0.01	0.02	0.03	0.02	0.01	0.01
23:14 00:14	0.00	0.00	-0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
00:14 01:14	0.00	-0.01	-0.01	0.01	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
01:14 02:14	0.00	-0.01	-0.01	0.00	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
02:14 03:14	0.00	-0.00	-0.01	0.00	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
03:14 04:14	0.00	0.00	-0.00	0.00	0.00	-0.00	0.01	0.02	0.02	0.01	0.00	-0.00
04:14 05:14	0.00	0.00	0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
05:14 06:14	0.00	-0.00	-0.00	0.00	0.00	-0.00	0.01	0.02	0.02	0.01	0.00	-0.00
06:14 07:14	0.00	-0.01	-0.01	0.00	0.01	-0.01	0.01	0.01	0.01	0.00	-0.00	-0.00
07:14 08:14	0.00	-0.01	-0.01	0.01	0.01	-0.01	0.01	0.01	0.01	0.00	-0.00	-0.00
08:14 09:14	0.00	0.00	-0.00	-0.00	0.01	-0.00	0.01	0.02	0.01	0.01	0.00	-0.00
Current day: 1 / Total 2 days												

ST Table

Time	I_ST mV	II_ST mV	III_ST mV	aVR_ST mV	aVL_ST mV	aVF_ST mV	V1_ST mV	V2_ST mV	V3_ST mV	V4_ST mV	V5_ST mV	V6_ST mV
09:14 10:14	0.01	0.01	0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
10:14 11:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.01
11:14 12:14	0.01	0.02	0.01	-0.01	0.00	0.01	0.01	0.03	0.03	0.02	0.01	0.01
12:14 13:14	0.01	0.02	0.01	-0.01	-0.00	0.01	0.01	0.03	0.02	0.02	0.01	0.01
13:14 14:14	0.00	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
14:14 15:14	0.00	0.01	0.01	-0.01	0.00	0.01	0.02	0.02	0.02	0.02	0.01	0.01
15:14 16:14	0.00	-0.00	-0.01	0.00	0.01	-0.01	0.01	0.02	0.02	0.01	0.00	-0.00
16:14 17:14	0.01	0.00	-0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
17:14 18:14	0.00	0.01	0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
18:14 19:14	0.01	0.01	0.01	-0.01	-0.00	0.01	0.01	0.03	0.02	0.02	0.01	0.01
19:14 20:14	0.00	0.01	0.01	-0.01	-0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
20:14 21:14	0.00	0.00	0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.00	0.00
21:14 22:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
22:14 23:14	0.00	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.01
23:14 00:14	0.01	0.02	0.01	-0.01	-0.00	0.01	0.01	0.02	0.02	0.02	0.01	0.01
00:14 01:14	0.00	0.01	0.00	-0.00	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
01:14 02:14	0.00	0.00	-0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
02:14 03:14	0.01	0.01	0.01	-0.01	0.00	0.01	0.01	0.02	0.02	0.02	0.01	0.01
03:14 04:14	0.00	-0.00	-0.00	0.00	0.01	-0.00	0.01	0.02	0.02	0.01	0.01	0.00
04:14 05:14	0.01	0.00	-0.00	-0.00	0.01	0.00	0.01	0.02	0.02	0.01	0.01	0.00
05:14 06:14	0.00	-0.01	-0.01	0.01	0.01	-0.01	0.01	0.02	0.01	0.01	-0.00	-0.01
06:14 07:14	0.00	0.00	-0.00	-0.00	0.00	-0.00	0.01	0.02	0.02	0.01	0.00	-0.00
07:14 08:14	0.01	0.01	0.00	-0.00	0.00	0.01	0.01	0.02	0.02	0.01	0.01	0.00
08:14 09:03	0.01	0.02	0.01	-0.01	-0.00	0.01	0.01	0.03	0.03	0.02	0.01	0.01
All	0.00	0.01	0.00	-0.00	0.00	0.00	0.01	0.02	0.02	0.01	0.01	0.00
Current day: 2 / Total 2 days												


HRV Chart


Time	Max R-R ms	Min R-R ms	Av. R-R ms	SDNN ms	HRVIndex ms	CV	MSD ms	rMSSD ms	PNN50 %	VLI ms	VAI degree	ULF ms*ms	VLF ms*ms	LF ms*ms	HF ms*ms
09:14 10:14	1035	330	660	73.6	156.8	0.09	22	52	8.2	89	0.92	1091.5	782.8	176.8	230.4
10:14 11:14	1195	435	712	80.4	95.3	0.11	32	71	13.2	91	1.30	422.7	1013.7	437.6	490.0
11:14 12:14	1185	455	727	72.4	94.4	0.10	31	72	10.2	73	1.16	132.5	556.2	468.4	636.5
12:14 13:14	1180	320	736	87.0	82.4	0.10	30	75	9.4	94	1.11	494.2	1037.6	633.9	752.1
13:14 14:14	1330	465	768	117.8	144.3	0.15	81	132	35.6	112	2.96	216.2	1134.7	1055.2	1656.5
14:14 15:14	1375	470	715	116.5	149.5	0.16	83	129	41.7	122	3.26	111.2	1209.4	1696.4	1745.3
15:14 16:14	1400	355	664	115.3	191.0	0.15	79	121	43.4	118	3.33	388.0	2017.1	1251.9	1426.5
16:14 17:14	1070	430	655	84.7	136.1	0.12	41	81	21.8	90	1.79	520.6	1098.4	455.7	606.6
17:14 18:14	1130	480	723	68.2	65.4	0.09	27	68	9.0	66	1.02	102.2	508.2	450.4	515.1
18:14 19:14	1105	365	717	52.6	63.2	0.07	17	49	4.1	55	0.67	78.7	394.7	244.7	312.5
19:14 20:14	1040	340	704	55.6	75.3	0.07	17	47	5.6	67	0.70	285.9	504.8	156.1	199.2
20:14 21:14	1125	470	740	70.7	90.6	0.09	34	78	11.1	62	1.24	97.6	470.4	314.1	450.3
21:14 22:14	1240	440	715	77.2	79.6	0.09	30	72	11.8	83	1.19	515.0	826.3	319.6	396.0
22:14 23:14	1210	470	718	86.4	71.2	0.12	63	105	31.5	70	2.44	37.6	340.4	782.0	1214.3
23:14 00:14	1285	305	727	113.2	110.6	0.15	92	134	47.9	104	3.67	94.3	935.8	1600.7	1645.6
00:14 01:14	1325	455	741	128.9	210.4	0.17	106	146	55.3	136	4.20	148.3	1522.9	2215.4	1961.4
01:14 02:14	1350	470	781	141.5	237.1	0.18	136	170	68.1	141	5.22	271.0	2702.9	2469.2	1737.6
02:14 03:14	1475	475	809	140.1	205.6	0.17	106	152	50.8	159	3.92	200.2	2156.0	2477.6	2040.9
03:14 04:14	1545	435	760	155.8	244.8	0.17	78	130	35.0	183	2.82	1340.0	4466.6	1684.2	1892.9
04:14 05:14	1400	465	819	133.0	223.2	0.16	95	145	42.1	154	3.35	186.2	2127.3	2286.8	1817.7
05:14 06:14	1465	465	795	141.6	243.4	0.18	88	137	42.9	167	3.21	283.8	2875.2	2588.2	1719.1
06:14 07:14	1410	405	679	129.4	205.7	0.13	44	89	21.4	160	1.78	2321.9	4028.8	1063.6	948.5
07:14 08:14	1095	420	654	79.8	108.9	0.11	33	71	15.4	91	1.43	445.2	1059.3	403.4	528.4
08:14 09:14	1170	325	694	71.0	98.3	0.07	19	47	6.6	92	0.82	774.2	838.9	203.9	220.7
Current day: 1 / Total 2 days															

HRV Chart

Time	Max R-R ms	Min R-R ms	Av. R-R ms	SDNN ms	HRVIndex ms	CV	MSD ms	rMSSD ms	PNN50 %	VLI ms	VAI degree	ULF ms*ms	VLf ms*ms	LF ms*ms	HF ms*ms
09:14 10:14	1080	340	686	67.1	88.8	0.08	20	48	7.5	84	0.85	603.0	764.1	196.1	219.7
10:14 11:14	1105	310	686	59.3	111.7	0.07	15	38	3.5	74	0.61	403.1	686.4	222.7	222.3
11:14 12:14	1160	520	734	56.3	75.4	0.07	17	47	3.1	61	0.64	101.8	440.0	354.9	352.6
12:14 13:14	1375	320	721	94.3	136.2	0.10	28	74	8.5	106	1.05	1084.9	1443.3	601.6	835.0
13:14 14:14	1370	345	710	83.3	124.2	0.10	24	62	7.1	98	0.92	648.7	1348.3	553.4	540.6
14:14 15:14	1310	455	749	112.8	172.4	0.15	71	123	31.2	110	2.66	151.5	1104.3	1145.7	1531.9
15:14 16:14	1240	420	650	104.1	117.6	0.14	70	111	38.6	102	3.02	620.2	1660.4	815.0	1088.6
16:14 17:14	1015	430	650	74.4	127.8	0.11	41	77	22.2	75	1.80	274.4	718.9	349.3	547.1
17:14 18:14	1155	420	686	86.5	89.7	0.12	45	88	20.2	85	1.84	375.9	978.1	318.9	765.5
18:14 19:14	1040	455	700	62.9	71.1	0.08	26	62	9.6	64	1.03	147.8	538.6	227.7	388.3
19:14 20:14	1095	435	695	63.9	101.3	0.08	19	48	6.0	80	0.80	535.8	665.5	172.3	229.7
20:14 21:14	1035	360	685	69.2	88.6	0.09	21	54	8.0	83	0.89	580.5	805.8	162.5	245.7
21:14 22:14	1100	300	696	70.2	85.9	0.09	24	60	8.4	78	0.94	503.1	718.6	223.9	403.5
22:14 23:14	1240	360	716	100.5	122.4	0.14	60	104	28.7	101	2.35	86.6	853.4	1274.0	1302.8
23:14 00:14	1220	435	706	97.2	104.9	0.13	64	107	31.2	88	2.49	99.9	647.4	680.9	1493.4
00:14 01:14	1370	435	760	123.2	196.3	0.16	66	117	30.2	136	2.41	198.3	1683.6	1921.3	1640.2
01:14 02:14	1330	335	767	133.6	204.3	0.17	73	126	32.2	149	2.65	303.9	2484.9	1802.6	1884.5
02:14 03:14	1305	455	762	117.6	143.0	0.15	106	149	50.6	108	4.05	148.9	1341.9	1658.8	1817.9
03:14 04:14	1420	445	758	141.1	249.4	0.18	134	169	65.8	135	5.19	382.7	3128.6	2397.5	1722.7
04:14 05:14	1400	410	783	142.9	215.0	0.17	134	170	64.0	140	5.04	438.0	3348.4	2494.1	1595.2
05:14 06:14	1415	415	697	156.9	218.4	0.17	107	149	52.6	160	4.20	2399.0	5846.3	1836.2	1278.9
06:14 07:14	1120	400	659	104.0	144.1	0.13	54	94	28.3	110	2.28	1134.7	1588.4	608.0	959.1
07:14 08:14	1385	445	732	111.8	179.2	0.13	51	100	21.7	127	1.91	752.5	2121.0	786.2	1079.5
08:14 09:03	1135	300	717	67.5	100.3	0.08	27	64	9.3	71	1.07	95.1	510.3	374.3	558.5
All	1545	300	715	105.7	161.3	0.12	50	97	22.6	117	1.96	736.8	2104.3	845.9	905.2
Current day: 2 / Total 2 days															

HRV Chart





Max R-R = 1545.0 ms

Min R-R = 300.0 ms

Av. R-R = 715.0 ms

SDNN = 105.7 ms


HRVIndex = 161.3 ms

TriangleIndex = 322.6 ms

CV = 0.12

MeanSD = 88.9 ms


SDANN = 57.2 ms




MSD = 50.2 ms

rMSSD = 97.1 ms

PNN50 = 22.6 %




Max R-R	1545.0 ms
Min R-R	300.0 ms
Av. R-R	715.0 ms
SDNN	105.7 ms
SDANN	57.2 ms
MeanSD	88.9 ms
CV	0.12
HRVIndex	161.3 ms
TriangleIndex	322.6 ms
rMSSD	97.1 ms
PNN50	22.6 %
VLI	116.8 ms
VAI	1.96 degree




Pa (-, +)	34.8 (%)
Pb (+, +)	22.8 (%)
Pc (-, -)	27.4 (%)
Pd (+, -)	15.0 (%)

N-N Frequence Domain Graph




ULFnoum = 16.0 %	VLFnorm = 45.8 %
LFnorm = 18.4 %	HFnorm = 19.7 %
LF/HF = 0.9	ULF = 736.8 ms*ms
VLF = 2104.3 ms*ms	LF = 845.9 ms*ms
HF = 905.2 ms*ms	total power = 4592.2 ms*ms


N-N Freq. Domain 3D Graph





Trends



Trends




Trends



V Run Table


ID	Start Time	Length	Av. HR	Beats
1	17:04:24	1690	142	4
2	20:19:04	1510	119	3
3	07:11:14	1345	133	3
4	06:13:22	1320	136	3
5	13:11:49	1290	139	3
6	16:27:52	1280	140	3
7	13:24:02	1280	140	3
8	11:01:45	1260	142	3
9	07:08:45	1255	143	3
10	11:01:00	1210	148	3
11	20:40:28	1175	153	3
12	07:03:44	1175	153	3
13	06:43:24	1090	165	3
14	06:43:18	1085	165	3
15	06:41:49	1080	166	3
Unit	---	ms	bpm	time
Current page: 1 / Total 1 pages				

DEMO REPORT 12 Lead ecg



2023-04-15 03:07:16				Clase: MinHR				HR: 63bpm				10 mm/mV 25.0 mm/s			
---------------------	--	--	--	--------------	--	--	--	-----------	--	--	--	--------------------	--	--	--

DEMO REPORT 12 Lead ecg



2023-04-15 03:38:43				Clase: MaxHR				HR: 121bpm				10 mm/mV 25.0 mm/s			
---------------------	--	--	--	--------------	--	--	--	------------	--	--	--	--------------------	--	--	--

Patient Name : **DEMO REPORT**

Gender : **MALE**

Age: **77**

Analysis Time
Duration: **6 Day 0 Hrs 3 Min**

Patient ID : **214943** Indication: **Dizziness**

Analysed Beats: **441195**

From: **2024/12/06@21:32**

Organization : **Arch Plus Healthcare Pvt L** Referred By: **Self**

To : **2024/12/12@21:40**

ATRIAL ECTOPICS

Total : 196|0.04%
Isolated : 155
Couplet : 3
Longest : 13 beats(133 bpm) 17:58:15 07-12
Fastest : 133 bpm (13 beats) 17:58:15 07-12
Trigeminy : 0%
Bigeminy : 0%

AFIB / FLUTTER

No AFIB Noted

VENTRICULAR ECTOPICS

Total : 184|0.04%
Isolated : 184
Couplet : 0
Longest :
Fastest :
Trigeminy : 0%
Bigeminy : 0.00%
Episodes : 0

PAUSES

Pauses in excess of 2 Sec : 0
Longest RR :


FINDINGS

- # Analysed an ECG recording of 6d:0 Hrs:3 Min:17 Sec.
- # Average HR was 61 bpm.
- # Minimum HR was 44 bpm at 22:54:00 | Maximum HR was 97 bpm at 10:46:26.
- # Ventricular ectopy was 184 (0.04 % of total analysed beats), with 0 V-pairs, 0 V-triplet 0 V-Run.
- # Ventricular Bigeminy events were 4 (0.00 % of total analysed beats) and Ventricular Trigeminy events were 0 (0 % of total analysed beats).
- # SupraVentricular ectopy was 196 (0.04 % of total analysed beats), with 3 SV-pairs, 4 SV-triplet 3 SV-Run. Longest SV run noted was of 13 beats on 2024-12-07 17:58:15.
- # Supraventricular run with maximum HR was 133 bpm (13 beats) on 2024-12-07 17:58:15.
- # No AV Block/ No Complete Heart Block/ No Ventricular Tachycardia/ No atrial fibrillation noted.
- # No Triggered Event Noted.



IMPRESSION: Sinus with First Degree AV Block, Intraventricular Conduction Delay (IVCD), PACs and PVCs were noted.

Atrial Couplets, Triplets, Runs and Ventricular Bigeminy noted.

Premature Ventricular Ectopics




Premature Atrial Ectopics



SR, 51.99 %
PVC, 0.04 %
PAC, 0.04 %
SB, 47.92 %


DAY 1 (2024-12-06 21:32 - 2024-12-07 21:32)

Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
21	1500	85	41	55	0	1488	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
22	3345	72	47	55	0	2981	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
23	3309	96	44	55	0	3060	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	300	65	53	56	0	276	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	3410	95	40	57	0	2990	0	3	0	0	0	0	0	0	0	0	0	0	0	0	0	0
2	3383	79	45	56	0	3184	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
3	3278	84	45	54	0	3132	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
4	3250	108	46	54	0	3108	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
5	3070	102	39	56	0	2780	0	3	0	0	0	0	2	0	0	0	0	0	0	0	0	0
6	3506	117	38	59	0	2557	0	2	0	0	0	0	4	0	1	0	0	0	0	0	0	0
7	3431	86	53	61	0	2076	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
8	3914	98	45	66	0	54	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	4076	130	48	70	0	295	0	4	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	3870	104	49	65	0	1138	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
11	4095	120	42	68	0	95	0	4	0	0	0	0	1	0	0	0	0	0	0	0	0	0
12	4093	120	42	68	0	169	0	3	0	0	0	0	4	0	1	0	0	0	0	0	0	0
13	4166	128	45	69	0	10	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
14	4063	118	46	67	0	14	0	3	0	0	0	0	1	0	0	0	0	0	0	0	0	0
15	3894	102	50	64	0	5	0	2	0	0	0	0	3	0	0	0	0	0	0	0	0	0
16	4554	118	50	77	0	3	0	0	0	0	0	0	6	0	0	0	0	0	0	0	0	0
17	4732	154	62	79	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0
18	4216	138	47	70	0	5	0	2	0	0	0	0	1	0	0	0	0	0	0	0	0	0
19	3856	92	51	64	0	178	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
20	3589	96	45	59	0	1485	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0




DAY 2 (2024-12-07 21:32 - 2024-12-08 21:32)

Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
22	1996	74	34	52	0	1978	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
23	3129	77	38	52	0	3009	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
0	2007	69	50	52	0	1980	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	3455	100	39	58	0	2005	0	1	0	0	0	0	3	0	0	0	0	0	0	0	1	0
2	3578	97	43	59	0	2043	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	3248	75	44	53	0	3194	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
4	3302	78	51	54	0	2996	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
5	3260	92	42	54	0	3152	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
6	3398	82	44	56	0	2899	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
7	3674	92	51	66	0	1092	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
8	3991	103	47	66	0	22	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	3968	100	49	67	0	275	0	2	0	0	0	0	1	0	0	0	0	0	0	0	0	0
10	4168	120	51	71	0	387	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
11	3622	96	41	60	0	2199	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
12	3753	106	44	62	0	1118	0	4	0	0	0	0	0	0	0	0	0	0	0	0	0	0
13	3799	103	40	63	0	883	0	5	0	0	0	0	1	0	0	0	0	0	0	0	0	0
14	3371	100	41	55	0	3137	0	2	0	0	0	0	1	0	0	0	0	0	0	0	0	0
15	3439	78	44	57	0	2700	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
16	3920	105	50	65	0	102	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
17	3770	119	51	62	0	201	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
18	3642	94	49	60	0	1367	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
19	3442	126	46	57	0	3140	0	1	0	0	0	0	2	1	0	0	0	0	0	0	0	0
20	3531	102	40	59	0	2318	0	1	0	0	0	0	0	1	0	0	0	0	0	0	0	0




DAY 3 (2024-12-08 21:32 - 2024-12-09 21:32)

Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
21	210	68	54	60	0	96	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
22	3367	103	41	56	0	2960	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
23	3154	166	45	52	0	2980	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	175	60	53	55	0	173	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	3148	82	43	52	0	3074	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
2	2481	71	50	55	0	2226	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	3334	74	43	55	0	3039	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
4	3318	106	37	55	0	3161	0	1	0	0	0	0	4	0	0	0	0	0	0	0	0	0
5	3396	83	45	56	0	2764	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
6	2477	117	49	59	0	1534	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0
7	3779	89	53	67	0	235	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
8	3985	125	51	69	0	121	0	2	0	0	0	0	4	0	0	0	0	0	0	0	0	0
9	4034	140	49	71	0	640	0	3	0	0	0	0	1	0	0	0	0	0	0	0	0	0
10	3694	122	42	62	0	1005	0	4	0	0	0	0	4	0	0	0	0	0	0	0	0	0
11	3811	135	46	63	0	381	0	2	0	0	0	0	1	0	0	0	0	0	0	0	0	0
12	3827	103	47	64	0	562	0	1	0	0	0	0	2	0	0	0	0	0	0	0	0	0
13	3789	104	45	63	0	599	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
14	3630	76	54	60	0	1360	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
15	3574	78	54	59	0	2222	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
16	3894	101	50	65	0	314	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
17	4408	135	54	74	0	3	0	2	0	0	0	0	2	0	0	0	0	0	0	0	0	0
18	4247	111	53	71	0	33	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
19	3704	91	50	62	0	330	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
20	3643	91	47	60	0	1228	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0




DAY 4 (2024-12-09 21:32 - 2024-12-10 21:32)

Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
22	2817	130	39	66	0	98	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
23	3524	131	40	59	0	2093	0	8	0	0	3	0	2	0	0	0	0	0	0	0	1	0
0	725	77	41	52	0	721	0	3	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	3236	82	44	54	0	2935	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
2	3064	72	45	53	0	2852	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	3274	95	43	54	0	2999	0	0	0	0	0	0	5	0	0	0	0	0	0	0	0	0
4	3524	75	54	59	0	2319	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
5	3242	81	43	55	0	2913	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
6	3080	85	50	55	0	2729	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
7	3500	81	48	58	0	2564	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
8	4032	122	47	67	0	17	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	559	89	53	64	0	12	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	3765	92	38	67	0	178	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
11	3742	101	43	62	0	871	0	1	0	0	0	0	2	0	0	0	0	0	0	0	0	0
12	3859	107	49	64	0	376	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
13	4098	120	53	68	0	9	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
14	4119	130	43	69	0	24	0	11	0	0	0	0	2	0	0	0	0	0	0	0	0	0
15	4109	127	44	69	0	11	0	13	0	0	0	0	0	0	0	0	0	0	0	0	0	0
16	3746	120	52	63	0	1106	0	5	0	0	0	0	1	0	0	0	0	0	0	0	0	0
17	4118	132	53	69	0	665	0	0	0	0	0	0	4	0	0	0	0	0	0	0	0	0
18	3725	110	42	62	0	818	0	8	0	0	0	0	2	0	0	0	0	0	0	0	0	0
19	3437	104	39	57	0	2807	0	3	0	0	0	0	3	0	0	0	0	0	0	0	0	0
20	3419	120	38	58	0	2501	0	1	0	0	0	0	3	1	0	0	0	0	0	0	0	0




DAY 5 (2024-12-10 21:32 - 2024-12-11 21:32)


Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
21	623	104	50	61	0	249	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
22	3205	106	38	53	0	3155	0	4	0	0	0	0	0	0	0	0	0	0	0	0	0	0
23	3202	78	43	53	0	3112	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	2207	61	50	51	0	2206	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	3157	70	48	52	0	3094	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
2	3085	60	44	51	0	3080	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	2956	71	48	53	0	2851	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
4	2815	78	51	54	0	2687	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
5	2434	94	43	53	0	2371	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
6	3402	88	48	58	0	2561	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0
7	3657	86	49	62	0	1293	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
8	3650	83	52	62	0	1539	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	3638	122	41	61	0	1808	0	4	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	4461	116	53	78	0	638	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
11	4460	107	50	75	0	2	0	0	0	0	0	0	3	0	0	0	0	0	0	0	0	0
12	3634	119	41	60	0	1885	0	3	0	0	0	0	0	0	0	0	0	0	0	0	0	0
13	3933	82	56	66	0	141	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
14	3614	104	43	60	0	1372	0	2	0	0	0	0	1	0	0	0	0	0	0	0	0	0
15	3563	84	52	60	0	1681	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
16	3457	108	48	57	0	2473	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
17	3590	93	43	60	0	1661	0	12	0	0	1	0	1	0	0	0	0	0	0	0	0	0
18	3494	73	49	58	0	2455	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
19	3409	79	51	57	0	2703	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
20	3462	120	40	57	0	2698	0	3	0	0	0	0	0	0	0	0	0	0	0	0	0	0



DAY 6 (2024-12-11 21:32 - 2024-12-12 21:32)

Hour	Beats	Max HR	Min HR	Avg HR	Tachy (SR > 100)	Brady (SR < 60)	SP > 2 sec	V ectopy	V Pair	V Triplet	V Bigemni	V Trigemny	A Ectopy	A Pair	A Triplet	A Bigemny	A Trigemny	VT	SVT	AFi b	Atrial Run	Ventricular Run
21	1404	96	44	56	0	177	0	1	0	0	0	0	1	0	0	0	0	0	0	0	0	0
22	3617	80	45	61	0	1996	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0
23	2305	100	45	65	0	175	0	3	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	2207	61	50	51	0	206	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	357	70	48	52	0	94	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
2	385	60	44	51	0	80	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	256	71	48	53	0	51	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
4	215	78	51	54	0	687	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
5	434	94	43	53	0	371	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
6	1402	88	48	58	0	261	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
7	367	86	49	62	0	123	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
8	650	83	52	62	0	139	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	338	122	41	61	0	108	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	1461	116	53	78	0	38	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
11	460	107	50	75	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
12	334	119	41	60	0	85	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
13	333	82	56	66	0	11	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
14	1614	104	43	60	0	72	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
15	363	84	52	60	0	181	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
16	347	108	48	57	0	73	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
17	1590	93	43	60	0	61	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
18	394	73	49	58	0	55	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
19	109	79	51	57	0	43	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
20	62	120	40	57	0	98	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0







ACTIVITY vs TIME

ACTIVITY vs ARRHYTHMIA




Lorenz Plot - All Beats


Lorenz Plot - Normal Beats





Maximum HR




Minimum HR




First degree AV block






IVCD














PVC Bigeminy




PVC









2024-12-09 21:48:12 1 mins




2024-12-09 22:49:30 1 mins




2024-12-09 23:31:42 1 mins




2024-12-10 14:52:44 1 mins




2024-12-10 15:03:27 1 mins




2024-12-10 15:15:15 1 mins




2024-12-10 15:15:42 1 mins




2024-12-10 15:37:46 1 mins




2024-12-10 15:39:13 1 mins







2024-12-10 16:06:06 1 mins




2024-12-10 22:04:34 1 mins






2024-12-11 09:00:28 1 mins











PAC Triplet






PAC Couplet




PAC




2024-12-07 16:29:43

1 mins




2024-12-07 16:29:57

1 mins






2024-12-07 18:12:24

1 mins





2024-12-07 21:54:53

1 mins




Sinus Bradycardia







2024-12-07 01:32:05 **1 mins**






2024-12-07 01:57:16 **1 mins**




2024-12-07 05:21:01 **1 mins**








2024-12-07 05:23:48 **1 mins**



Normal Sinus Rhythm





ARCHPLUS HEALTHCARE PVT. LTD.

Patient Name:	DEMO REPORT	www.aphcpl.com	+91-9723111128	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A		
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M		



Patient Information

Address:	VASTRAPUR	Height:	165cm
		Weight:	64kg
Out Patient No.:	N/A	Nationality:	INDIAN
Admission No.:	N/A	Telephone:	N/A
Email:			

Current Medications

Examine Result

All BP Averages:	127/63mmHg	BP threshold:	130/80mmHg
Day BP Averages:	123/62mmHg	BP threshold:	135/85mmHg
Night BP Averages:	137/67mmHg	BP threshold:	120/70mmHg

Day BP Load Value:Nomal<40%	Night BP Load Value:Nomal<50%
SYS(>135mmHg) 33.6%	SYS(>120mmHg) 85.7%
DIA(>85mmHg) 3.4%	DIA(>70mmHg) 37.1%

Maximum SYS	178mmHg	Time	2024/1/1 08:45	Minimum SYS	66mmHg	Time	2023/12/30 11:15
Maximum DIA	99mmHg	Time	2024/1/1 08:45	Minimum DIA	30mmHg	Time	2023/12/30 19:15

Circadian rhythm of BP:SYS Night Dec.	-11.3%	DIA Night Dec.	-8.9%	Normal:10%-20%
--	--------	-----------------------	-------	-----------------------

BP CV:	All:SYS	18.3%	DIA	21.0%	Atherogenic Index:	0.5
	Day:SYS	19.4%	DIA	22.2%	24 Hours SI:	SYS 8.7 DIA 7.8
	Night:SYS	11.8%	DIA	15.3%	Morning Surge:	15 mmHg

Comments & Diagnoses

ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Decision Statistics

Samples used/attempted: 154/189 (81.5%)

BP Statistics

Overall Statistics			Total Samples Used:		154			
	Maximum	Time	Minimum	Time	Averages	SD	SE	CV
Systolic(mmHg)	178	08:45	66	11:15	127	23.2	1.9	0.2
Diatolic(mmHg)	99	08:45	30	19:15	63	13.2	1.1	0.2
Pulse Rate(BPM)	85	13:48	58	12:00	70	5.0	0.4	<0.1
MAP(mmHg)	132	16:45	42	11:15	87	18.1	1.5	0.2
PP(mmHg)	91	23:03	24	10:46	64	--	--	--
Systolic> 130mmHg		50.6%		Diatolic> 80mmHg		7.1%		
Awake Statistics			Total Samples Used:		119			
	Maximum	Time	Minimum	Time	Averages	SD	SE	CV
Systolic(mmHg)	178	08:45	66	11:15	123	23.9	2.2	0.2
Diatolic(mmHg)	99	08:45	30	19:15	62	13.7	1.3	0.2
Pulse Rate(BPM)	85	13:48	58	12:00	70	5.2	0.5	<0.1
MAP(mmHg)	132	16:45	42	11:15	84	18.1	1.7	0.2
PP(mmHg)	86	07:15	24	10:46	61	--	--	--
Systolic> 135mmHg		33.6%		Diatolic> 85mmHg		3.4%		
Sleep Statistics			Total Samples Used:		35			
	Maximum	Time	Minimum	Time	Averages	SD	SE	CV
Systolic(mmHg)	161	23:03	96	04:00	137	16.2	2.7	0.1
Diatolic(mmHg)	90	05:00	39	23:00	67	10.3	1.7	0.2
Pulse Rate(BPM)	77	22:30	59	03:30	69	4.2	0.7	<0.1
MAP(mmHg)	131	05:00	67	23:00	98	13.4	2.3	0.1
PP(mmHg)	91	23:03	44	04:00	70	--	--	--
Systolic> 120mmHg		85.7%		Diatolic> 70mmHg		37.1%		

ARCHPLUS HEALTHCARE PVT. LTD. www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Ambulatory Blood Pressure Summary

Length of study:	47H57M
Number of successful readings:	154

Blood Pressure Load

Systolic BP threshold awake:	135 mmHg	Above threshold: 33.6%
Diatolic BP threshold awake:	85 mmHg	Above threshold: 3.4%
Systolic BP threshold asleep:	120 mmHg	Above threshold: 85.7%
Diatolic BP threshold asleep:	70 mmHg	Above threshold: 37.1%

Averages

Awake BP	123/62mmHg	Awake PP	61 mmHg	Awake pulse rate	70BPM
Asleep BP	137/67mmHg	Asleep PP	70 mmHg	Asleep pulse rate	69BPM


Clinical Interpretation

ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

BP Trend vs Time




ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

BP Trend vs Time




ARCHPLUS HEALTHCARE PVT. LTD.


www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M


Histogram(All)



Systolic(mmHg)



Diastolic(mmHg)



Pulse Rate(BPM)

ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128



Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Pie Chart(All)

Maximum: 178(mmHg)(08:45)

Minimum: 66(mmHg)(11:15)

Averages: 127(mmHg)





Systolic(mmHg)

Maximum: 99(mmHg)(08:45)

Minimum: 30(mmHg)(19:15)

Averages: 63(mmHg)





Diastolic(mmHg)

Maximum: 85(BPM)(13:48)

Minimum: 58(BPM)(12:00)

Averages: 70(BPM)




Pulse Rate(BPM)

ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Correlation line







ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M




Frequency Histograms



ARCHPLUS HEALTHCARE PVT. LTD.

www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M



ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Measured Data									
Num	Date	Time	SYS	MAP	DIA	PP	PR	State	Comment
1+	2023/12/30	10:05	117	90	58	59	63	0	
2	2023/12/30	10:15	96	67	55	41	63	0	
3	2023/12/30	10:30	141 ↑	96	69	72	69	0	
4	2023/12/30	10:45	116	72	51	65	66	0	
5	2023/12/30	11:00	90	63	47	43	62	0	
6	2023/12/30	11:15	66	42	31	35	62	0	
7	2023/12/30	11:31	71	50	39	32	61	0	
8	2023/12/30	11:45	97	63	47	50	60	0	
9	2023/12/30	12:00	107	72	50	57	58	0	
10	2023/12/30	12:15	107	72	56	51	59	0	
11	2023/12/30	12:30	146 ↑	94	73	73	66	0	
12	2023/12/30	12:45	139 ↑	108	66	73	66	0	
13	2023/12/30	13:00	131	102	78	53	68	0	
14	2023/12/30	13:15	133	82	58	75	67	0	
15	2023/12/30	13:30	135	81	62	73	68	0	
16	2023/12/30	13:45	160 ↑	103	81	79	68	0	
17	2023/12/30	14:00	153 ↑	94	73	80	69	0	
18	2023/12/30	14:15	134	101	70	64	66	0	
19++	2023/12/30	14:37	87	55	46	41	66	0	
20	2023/12/30	14:46	100	71	55	45	68	0	
21	2023/12/30	15:00	92	64	49	43	68	0	
22++	2023/12/30	15:23	135	85	69	66	76	0	
23	2023/12/30	15:30	148 ↑	95	72	76	73	0	
24	2023/12/30	15:45	109	60	42	67	66	0	
25	2023/12/30	16:00	146 ↑	90	68	78	71	0	
26	2023/12/30	16:16	143 ↑	105	77	66	74	0	
27	2023/12/30	16:30	159 ↑	101	80	79	74	0	
28++	2023/12/30	16:51	113	77	58	55	66	0	
29	2023/12/30	17:00	103	73	54	49	68	0	
30	2023/12/30	17:15	95	63	48	47	66	0	
31++	2023/12/30	17:33	101	71	54	47	68	0	
32	2023/12/30	17:45	86	56	31	55	68	0	
33	2023/12/30	18:01	84	58	42	42	68	0	
34	2023/12/30	18:15	86	60	45	41	70	0	
35	2023/12/30	18:30	95	69	49	46	69	0	
36	2023/12/30	18:45	94	67	49	45	68	0	
37	2023/12/30	19:00	90	57	43	47	67	0	
38	2023/12/30	19:15	97	46	30	67	68	0	
39	2023/12/30	19:30	104	74	54	50	68	0	
40	2023/12/30	19:46	109	71	54	55	70	0	

ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Measured Data									
Num	Date	Time	SYS	MAP	DIA	PP	PR	State	Comment
42	2023/12/30	20:15	91	63	48	43	68	0	
43++	2023/12/30	20:33	100	64	49	51	67	0	
44++	2023/12/30	20:48	117	92	82	35	68	0	
45	2023/12/30	21:00	130	85	64	66	71	0	
46	2023/12/30	21:15	131	87	59	72	69	0	
47++	2023/12/30	21:33	82	49	35	47	70	0	
48	2023/12/30	22:00	141 ↑	99	76	65	85	0	
49	2023/12/30	22:30	134 ↑	94	63	71	77	0	
50	2023/12/30	23:00	108	67	39	69	70	0	
51	2023/12/30	23:30	131 ↑	97	61	70	72	0	
52	2023/12/31	00:00	135 ↑	110	67	68	74	0	
53	2023/12/31	00:30	139 ↑	96	72 ↑	67	75	0	
54	2023/12/31	01:00	123 ↑	94	69	54	74	0	
55	2023/12/31	01:30	129 ↑	93	62	67	71	0	
56	2023/12/31	02:00	144 ↑	97	71 ↑	73	73	0	
57	2023/12/31	02:30	131 ↑	82	57	74	68	0	
58	2023/12/31	03:00	140 ↑	99	69	71	73	0	
59	2023/12/31	03:30	138 ↑	102	74 ↑	64	67	0	
60	2023/12/31	04:00	96	71	52	44	61	0	
61	2023/12/31	04:30	131 ↑	100	64	67	63	0	
62	2023/12/31	05:00	151 ↑	92	69	82	74	0	
63	2023/12/31	05:30	116	90	62	54	64	0	
64	2023/12/31	06:00	146 ↑	92	69	77	67	0	
65	2023/12/31	06:30	151 ↑	92	70	81	74	0	
66	2023/12/31	07:00	145 ↑	81	67	78	67	0	
67	2023/12/31	07:15	132	82	56	76	64	0	
68++	2023/12/31	07:33	128	77	59	69	66	0	
69	2023/12/31	08:01	111	85	61	50	70	0	
70	2023/12/31	08:15	118	88	60	58	72	0	
71	2023/12/31	08:30	106	76	56	50	71	0	
72	2023/12/31	08:45	107	73	58	49	69	0	
73	2023/12/31	09:00	117	79	63	54	68	0	
74++	2023/12/31	09:20	127	95	63	64	69	0	
75	2023/12/31	09:30	149 ↑	106	64	85	69	0	
76	2023/12/31	09:45	132	89	68	64	68	0	
77	2023/12/31	10:00	135	90	66	69	67	0	
78++	2023/12/31	10:18	130	88	74	56	68	0	
79	2023/12/31	10:30	118	78	56	62	68	0	
80+	2023/12/31	10:46	105	86	81	24	69	0	
81	2023/12/31	11:00	78	58	45	33	68	0	

ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Measured Data									
Num	Date	Time	SYS	MAP	DIA	PP	PR	State	Comment
83	2023/12/31	11:30	127	67	49	78	75	0	
84	2023/12/31	11:45	135	82	66	69	76	0	
85++	2023/12/31	12:03	149 ↑	102	76	73	77	0	
86	2023/12/31	12:15	149 ↑	83	69	80	78	0	
87	2023/12/31	12:30	170 ↑	113	86 ↑	84	83	0	
88	2023/12/31	12:45	126	95	76	50	82	0	
89	2023/12/31	13:00	128	114	99 ↑	29	78	0	
90	2023/12/31	13:15	156 ↑	98	76	80	78	0	
91	2023/12/31	13:30	151 ↑	100	66	85	78	0	
92++	2023/12/31	13:48	166 ↑	113	87 ↑	79	85	0	
93	2023/12/31	14:00	149 ↑	94	70	79	82	0	
94	2023/12/31	14:16	149 ↑	120	79	70	80	0	
95++	2023/12/31	14:33	138 ↑	93	74	64	80	0	
96	2023/12/31	14:45	87	61	47	40	74	0	
97	2023/12/31	15:00	96	59	46	50	72	0	
98	2023/12/31	15:15	106	70	54	52	72	0	
99	2023/12/31	15:30	136 ↑	84	67	69	73	0	
100	2023/12/31	15:45	140 ↑	93	73	67	77	0	
101	2023/12/31	16:00	148 ↑	104	69	79	79	0	
102	2023/12/31	16:15	154 ↑	97	73	81	79	0	
103	2023/12/31	16:30	137 ↑	102	70	67	73	0	
104	2023/12/31	16:45	159 ↑	132	74	85	78	0	
105	2023/12/31	17:00	138 ↑	89	66	72	73	0	
106	2023/12/31	17:15	134	101	74	60	75	0	
107	2023/12/31	17:45	93	55	40	53	70	0	
108	2023/12/31	18:00	108	72	56	52	70	0	
109	2023/12/31	18:16	101	73	48	53	70	0	
110	2023/12/31	18:30	97	64	42	55	69	0	
111	2023/12/31	18:45	121	75	58	63	69	0	
112	2023/12/31	19:00	117	80	63	54	68	0	
113	2023/12/31	19:15	126	79	59	67	68	0	
114	2023/12/31	19:30	140 ↑	82	65	75	65	0	
115	2023/12/31	19:45	113	82	45	68	68	0	
116	2023/12/31	20:00	126	81	48	78	75	0	
117++	2023/12/31	20:18	135	89	65	70	70	0	
118	2023/12/31	20:31	120	97	79	41	67	0	
119	2023/12/31	20:45	126	76	59	67	69	0	
120	2023/12/31	21:00	132	101	57	75	71	0	
121	2023/12/31	21:16	108	77	53	55	75	0	
122	2023/12/31	21:30	115	89	57	58	67	0	

ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Measured Data									
Num	Date	Time	SYS	MAP	DIA	PP	PR	State	Comment
124	2023/12/31	22:01	130 ↑	102	61	69	67	0	
125	2023/12/31	22:30	153 ↑	100	79 ↑	74	68	0	
126++	2023/12/31	23:03	161 ↑	122	70	91	67	0	
127	2023/12/31	23:30	145 ↑	99	74 ↑	71	68	0	
128++	2024/1/1	00:03	158 ↑	103	83 ↑	75	68	0	
129	2024/1/1	00:30	144 ↑	99	79 ↑	65	65	0	
130++	2024/1/1	01:07	130 ↑	104	61	69	72	0	
131	2024/1/1	01:31	151 ↑	115	80 ↑	71	66	0	
132	2024/1/1	02:00	127 ↑	102	62	65	62	0	
133	2024/1/1	02:30	102	72	52	50	64	0	
134	2024/1/1	03:01	133 ↑	95	47	86	66	0	
135	2024/1/1	03:30	115	89	57	58	59	0	
136++	2024/1/1	04:06	137 ↑	90	66	71	70	0	
137	2024/1/1	04:31	156 ↑	102	79 ↑	77	73	0	
138	2024/1/1	05:00	156 ↑	131	90 ↑	66	68	0	
139	2024/1/1	05:30	157 ↑	91	72 ↑	85	68	0	
140	2024/1/1	06:00	159 ↑	127	73 ↑	86	70	0	
141	2024/1/1	06:30	150 ↑	109	72 ↑	78	71	0	
142	2024/1/1	07:00	155 ↑	125	77	78	68	0	
143	2024/1/1	07:15	162 ↑	126	76	86	68	0	
144	2024/1/1	07:31	130	91	70	60	71	0	
145	2024/1/1	07:45	110	80	64	46	70	0	
146	2024/1/1	08:00	152 ↑	90	70	82	74	0	
147	2024/1/1	08:15	153 ↑	106	84	69	76	0	
148	2024/1/1	08:30	154 ↑	96	75	79	72	0	
149	2024/1/1	08:45	178 ↑	119	99 ↑	79	75	0	
150	2024/1/1	09:00	166 ↑	103	82	84	72	0	
151++	2024/1/1	09:19	144 ↑	109	70	74	68	0	
152	2024/1/1	09:31	139 ↑	102	58	81	66	0	
153	2024/1/1	09:45	140 ↑	79	65	75	64	0	
154	2024/1/1	10:02	137 ↑	83	66	71	70	0	

ARCHPLUS HEALTHCARE PVT. LTD.					
www.aphcpl.com +91-9723111128					
Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Hourly Summary Report

Time	SYS	DIA	MAP	PR	RPP	Number	Duration
10:05(1)	118	58	81	65	7670	4	55
11:00(1)	81	41	55	61	4941	4	60
12:00(1)	125	61	87	62	7750	4	60
13:00(1)	140	70	92	68	9520	4	60
14:00(1)	119	61	80	67	7973	4	60
15:00(1)	121	58	76	71	8591	4	60
16:00(1)	140	71	93	71	9940	4	60
17:00(1)	96	47	66	68	6528	4	60
18:00(1)	90	46	64	69	6210	4	60
19:00(1)	100	45	62	68	6800	4	60
20:00(1)	101	57	72	68	6868	4	60
21:00(1)	114	53	74	70	7980	3	60
22:00(1)	138	70	97	81	11178	2	60
23:00(1)	120	50	82	71	8520	2	60
00:00(2)	137	70	103	75	10275	2	60
01:00(2)	126	66	94	73	9198	2	60
02:00(2)	138	64	90	71	9798	2	60
03:00(2)	139	72	101	70	9730	2	60
04:00(2)	114	58	86	62	7068	2	60
05:00(2)	134	66	91	69	9246	2	60
06:00(2)	149	70	92	71	10579	2	60
07:00(2)	135	61	80	66	8910	3	60
08:00(2)	111	59	81	71	7881	4	60
09:00(2)	131	65	92	69	9039	4	60
10:00(2)	122	69	86	68	8296	4	60
11:00(2)	113	54	73	73	8249	4	60
12:00(2)	149	77	98	80	11920	4	60
13:00(2)	150	82	106	80	12000	4	60
14:00(2)	131	68	92	79	10349	4	60
15:00(2)	120	60	77	74	8880	4	60
16:00(2)	150	72	109	77	11550	4	60
17:00(2)	122	60	82	73	8906	3	60
18:00(2)	107	51	71	70	7490	4	60
19:00(2)	124	58	81	67	8308	4	60
20:00(2)	127	63	86	70	8890	4	60
21:00(2)	116	54	83	70	8120	4	60
22:00(2)	142	70	101	68	9656	2	60

ARCHPLUS HEALTHCARE PVT. LTD.www.aphcpl.com +91-9723111128

Patient Name:	DEMO REPORT	Male/Female:	Male	Age:	81
Patient ID:	2312291816	Bed No.:	AT HOME	Department No.:	N/A
Test Begin:	2023/12/30 10:05	Test End:	2024/01/01 10:02	Duration:	47H57M

Hourly Summary Report

Time	SYS	DIA	MAP	PR	RPP	Number	Duration
23:00(2)	153	72	111	68	10404	2	60
00:00(3)	151	81	101	67	10117	2	60
01:00(3)	141	71	110	69	9729	2	60
02:00(3)	115	57	87	63	7245	2	60
03:00(3)	124	52	92	63	7812	2	60
04:00(3)	147	73	96	72	10584	2	60
05:00(3)	157	81	111	68	10676	2	60
06:00(3)	155	73	118	71	11005	2	60
07:00(3)	139	72	106	69	9591	4	60
08:00(3)	159	82	103	74	11766	4	60
09:00(3)	145	68	95	68	9860	5	62

PORTABLE NCV REORT

Full Name:	DEMO REPORT	Patient ID:	
Age:	46 year	Gender:	Male
Referring Dr. :		Reporting Dr.	
Technician :		Test Date:	23-06-2026 11:37
Diagnoses:			

NCS REPORT

IMPRESSION

MNCV

- Normal distal motor latencies, normal CMAPs amplitudes and normal conduction velocities in bilateral Median (APB), Ulnar (ADM), Peroneal (EDB) & bilateral Tibial (AH) nerves.

SNCV

- Normal sensory latencies, normal SNAPs amplitudes and normal sensory conduction velocities in right Median (Dig II), Ulnar (Dig V) & bilateral Sural nerves.
- Normal sensory latency, normal SNAPs amplitudes and mild reduced sensory conduction velocity in left Median (Dig II) nerve. nerves.

Late Response

- Normal F-waves latencies all tested nerves.

CONCLUSION

Above NCS study of all four limbs is does not show any significant abnormalities.
Please correlate clinically.

Dr.

MD, DNB Neurology
Consultant, Neurologist.

Motor Nerve Conduction









Nerve / Sites	Latency ms	Amp mV	Area mVms	Dur. ms	% Amp %	Segments	Dist. mm	d Lat. ms	Vel. m/s
R Median - APB									
Wrist	3.1	15.9	57.1	6.6	100	Wrist - APB			
Elbow	7.8	14.9	56.7	7.0	93.5	Elbow - Wrist	250	4.7	53
L Median - APB									
Wrist	3.8	10.5	45.4	7.1	100	Wrist - APB			
Elbow	8.5	10.3	45.4	7.0	97.4	Elbow - Wrist	250	4.8	52
R Ulnar - ADM									
Wrist	2.5	9.2	30.1	6.2	100	ADM - Wrist			
B.Elbow	7.8	8.7	28.5	6.2	94.9	B.Elbow - Wrist	260	5.2	50
L Ulnar - ADM									
Wrist	2.5	9.7	32.6	6.6	100	ADM - Wrist			
B.Elbow	7.4	9.3	31.0	7.0	95.4	B.Elbow - Wrist	260	4.9	54
R Peroneal - EDB									
Ankle	3.1	5.9	18.1	5.6	100	Ankle - EDB			
Fib head	11.8	4.7	15.4	5.8	79.1	Fib head - Ankle	380	8.7	44
L Peroneal - EDB									
Ankle	3.3	5.1	19.6	7.2	100	Ankle - EDB			
Fib head	12.8	3.8	14.9	6.6	74.4	Fib head - Ankle	380	9.5	40
R Tibial - AH									
Ankle	3.9	12.2	33.6	6.6	100	Ankle - AH			
Pop fossa	14.6	7.0	21.3	7.1	57.8	Pop fossa - Ankle	440	10.7	41
L Tibial - AH									
Ankle	4.7	8.7	30.3	6.6	100	Ankle - AH			
Pop fossa	15.0	5.0	14.5	6.9	57.6	Pop fossa - Ankle	440	10.3	43


Sensory Nerve Conduction


Nerve / Sites	Rec. Site	Onset Lat ms	P-P Amp μ V	Segments	Distance mm	d Lat. ms	Velocity m/s
R Median - Digit II (Antidromic)							
Wrist	Dig II	2.6	28.0	Wrist - Dig II	130	2.6	50
L Median - Digit II (Antidromic)							
Wrist	Dig II	3.2	38.2	Wrist - Dig II	130	3.2	41
R Ulnar - Digit V (Antidromic)							
Wrist	Dig V	2.3	24.3	Wrist - Dig V	110	2.3	48
L Ulnar - Digit V (Antidromic)							
Wrist	Dig V	2.3	36.3	Wrist - Dig V	110	2.3	49
R Sural - Ankle (Calf)							
1	Ankle	2.4	14.4	1 - Ankle	100	2.4	41
L Sural - Ankle (Calf)							
1	Ankle	2.5	13.8	1 - Ankle	100	2.5	40

F-Wave

Nerve	Fmin ms	FMax ms
R Median - APB	28.0	30.5
R Ulnar - ADM	28.7	30.9
R Peroneal - EDB	51.8	56.1
R Tibial - AH	55.7	58.8
L Median - APB	30.2	31.4
L Ulnar - ADM	29.6	31.9
L Peroneal - EDB	51.2	55.9
L Tibial - AH	56.7	58.8








PORTABLE VEP REPORT

Full Name:	DEMO REPORT	Patient ID:	
Age:	43 Years	Gender:	Female
Referring Dr. :		Reporting Dr.	
Technician :		Test Date:	28-11-2025 12:49
Diagnoses:			

VEP(Pattern) REPORT

INTERPRETATION:



Run	Side	Label	N75 ms	P100 ms	N145 ms	Ampl µV	Notes
2	Left	Oz-Fz	87.5	142.1	210.0	10.3	
1	Right	Oz-Fz	84.6	121.7	190.0	11.5	

Left eye showed P-100 latency of 142.1 ms with normal amplitude.

Right eye showed P-100 latency of 121.7 ms with normal amplitude.

CONCLUSION:

Above VEP study is shows prolonged P-100 latencies on bilateral Eyes with normal amplitude. This study is suggestive of bilateral anterior visual pathway dysfunction on both side. (Left >> Right)


Dr.

MD, DNB Neurology

VEP-Pattern - 1 Ch

Run	Side	Label	N75 ms	P100 ms	N145 ms	Ampl μV	Notes
Pattern VEP							
1	Left	Oz-Cz	87.5	142.1	210.0	10.3	
2	Right	Oz-Cz	84.6	121.7	190.0	11.5	

Run			▲ N75 ms	▲ P100 ms	▲ N145 ms	▲ Ampl μV
Pattern VEP						
1-2			2.9	20.4	20.0	1.1



Glucose Monitoring Report

Basic info

Name: DEMO REPORT

Type of Diabetes: Other

Device SN Code: 88F016FCD604A

Gender: Unknown

Age: --

Device Model: i3

Glucose Overview

Reporting period: 15/09/2026 ~ 23/09/2026


Monitoring Days: 8 Days

Time CGM Active: 100%

Glucose Metrics


Mean Glucose (MG)¹
Goal: <154mg/dL **113 mg/dL**Glucose Management Indicator (GMI)²
Goal: <7% **6 %**Coefficient of variation (CV)³
Goal: ≤36% **23.7 %**

Time in Ranges⁴



Each 1% time in range = about 15 minutes


Ambulatory Glucose Profile (AGP)



Low Glucose Condition

Hypoglycemia event(<70mg/dL)⁵ **6** times Average Duration **59** minute(s)


Serious hypoglycaemia event(<54mg/dL)⁶ **2** times Average Duration **47** minute(s)




High Glucose Condition

Hyperglycaemia Event (>180mg/dL)⁷ **1** time(s) Average Duration **93** minute(s)

Serious Hyperglycaemia Event (>250mg/dL)⁸ **0** times



Multi-day Glucose Details



Date	15/09	16/09	17/09	18/09	19/09	20/09	21/09	22/09	23/09							
Number of glucose value	68	480	480	480	480	480	480	480	291							
Highest mg/dL	266	180	169	176	150	150	158	178	150							
Lowest mg/dL	144	75	55	48	71	42	57	78	80							
MG mg/dL	180	124	122	115	103	102	105	118	107							
TIR	61.8%	100%	97.5%	95.8%	100%	95%	92.9%	100%	100%							
TAR	38.2%	0%	0%	0%	0%	0%	0%	0%	0%							
TBR	0%	0%	2.5%	4.2%	0%	5%	7.1%	0%	0%							
LAGE ¹² mg/dL	122	105	114	128	79	108	101	100	70							
MAGE ¹³ mg/dL	--*	34	40	44	23	42	46	40	--*							
MODD ¹⁴ mg/dL	--*	--*	19	13	25	18	12	24	--*							
SDBG ¹⁵ mg/dL	28	24	25	28	20	21	23	25	24							
CV	15.6%	19.6%	20.3%	24.3%	19.3%	20.9%	21.7%	21.2%	22.3%							

* The data is less than 24 h and cannot be calculated.

| Multi-day segmented glucose details


Time period	00:00 ~ 01:00	01:00 ~ 02:00	02:00 ~ 03:00	03:00 ~ 04:00	04:00 ~ 05:00	05:00 ~ 06:00	06:00 ~ 07:00	07:00 ~ 08:00	08:00 ~ 09:00	09:00 ~ 10:00	10:00 ~ 11:00	11:00 ~ 12:00
Number of glucose value	160	160	160	160	160	160	160	160	160	160	160	160
Highest mg/dL	148	147	139	106	106	106	119	162	175	180	156	163
Lowest mg/dL	86	77	72	71	74	74	82	91	113	109	112	109
TIR	100%	100%	100%	100%	100%	100%	100%	100%	100%	100%	100%	100%
TAR	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%
TBR	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%	0%
MG mg/dL	113	110	98	91	90	91	96	108	148	142	132	133
CV	14.1%	18.2%	14.7%	9.8%	8.5%	6.9%	7.2%	12.5%	9.6%	11.8%	8.6%	11.7%
SDBG ¹⁵ mg/dL	16	20	14	9	8	6	7	14	14	17	11	16
Median line ¹⁶ mg/dL	110	105	99	92	90	92	96	110	137	139	132	127
IQR ¹⁷ 25%-75% range mg/dL	23	26	21	13	7	7	8	15	23	21	19	25
IVR ¹⁸ 5%-95% range mg/dL	58	58	43	30	26	22	25	38	43	52	44	46

Time period	12:00 ~ 13:00	13:00 ~ 14:00	14:00 ~ 15:00	15:00 ~ 16:00	16:00 ~ 17:00	17:00 ~ 18:00	18:00 ~ 19:00	19:00 ~ 20:00	20:00 ~ 21:00	21:00 ~ 22:00	22:00 ~ 23:00	23:00 ~ 00:00
Number of glucose value	160	160	151	140	140	140	140	140	148	160	160	160
Highest mg/dL	161	128	130	157	166	164	176	178	266	215	177	169
Lowest mg/dL	81	42	44	68	89	102	88	83	75	67	75	81
TIR	100%	74.3%	78.2%	94.3%	100%	100%	100%	100%	94.6%	83.7%	100%	100%
TAR	0%	0%	0%	0%	0%	0%	0%	0%	5.4%	11.3%	0%	0%
TBR	0%	25.7%	21.8%	5.7%	0%	0%	0%	0%	0%	5%	0%	0%
MG mg/dL	116	86	85	100	122	131	134	125	121	120	119	116
CV	13.7%	25.1%	23.8%	24.2%	15.9%	15.1%	17.8%	22.3%	30.7%	29.3%	23.6%	17.1%
SDBG ¹⁵ mg/dL	16	21	20	24	19	20	24	28	37	35	28	20
Median line ¹⁶ mg/dL	112	89	87	96	116	131	129	116	110	120	123	115
IQR ¹⁷ 25%-75% range mg/dL	25	33	33	31	33	38	40	47	45	46	40	26
IVR ¹⁸ 5%-95% range mg/dL	50	61	68	76	62	58	73	86	114	117	96	70


Daily Glucose Profiles

Fingerstick blood Exercise Diet Medication Insulin Multi Log Other


TAR TIR TBR




MG	Highest	Lowest
180mg/dL	266mg/dL	144mg/dL
TAR	TIR	TBR
38.2%	61.8%	--%
1h18min	2h6min	0min




MG	Highest	Lowest
124mg/dL	180mg/dL	75mg/dL
TAR	TIR	TBR
--%	100%	--%
0min	24h	0min




MG	Highest	Lowest
122mg/dL	169mg/dL	55mg/dL
TAR	TIR	TBR
--%	97.5%	2.5%
0min	23h24min	36min




MG	Highest	Lowest
115mg/dL	176mg/dL	48mg/dL
TAR	TIR	TBR
--%	95.8%	4.2%
0min	23h	1h




MG	Highest	Lowest
103mg/dL	150mg/dL	71mg/dL
TAR	TIR	TBR
--%	100%	--%
0min	24h	0min




MG	Highest	Lowest
102mg/dL	150mg/dL	42mg/dL
TAR	TIR	TBR
--%	95%	5%
0min	22h48min	1h12min



MG	Highest	Lowest
105 mg/dL	158 mg/dL	57 mg/dL
TAR	TIR	TBR
--%	92.9 %	7.1 %
0min	22h18min	1h42min



MG	Highest	Lowest
118 mg/dL	178 mg/dL	78 mg/dL
TAR	TIR	TBR
--%	100 %	--%
0min	24h	0min



MG	Highest	Lowest
107 mg/dL	150 mg/dL	80 mg/dL
TAR	TIR	TBR
--%	100 %	--%
0min	14h33min	0min


Explanation of metrics

1. Mean Glucose (MG): Refers to the average of all glucose measurements during the monitoring period to evaluate the overall glucose level.
2. Glucose Management Indicator (GMI): GMI is derived from at least 15 days of CGM data. The GMI may be similar to, higher than, or lower than the laboratory HbA1c.
3. Glucose Variability (CV): The degree of dispersion of all glucose measurements during the monitoring period is more suitable for comparing glucose variability in different populations and for evaluating glucose fluctuations.
4. Time in Ranges (TIR,TAR,TBR): Refers to the percentage of time during the monitoring period when glucose was within the target range, above the target range (distinguishing between high glucose and very high glucose), and below the target range (distinguishing between low glucose and very low glucose), respectively, and is used to evaluate the level of glycaemic control.
5. Hypoglycemia Event: Event starts when the glucose value is lower than the low limit of the target range (e.g., 3.9 mmol/L or 70 mg/dL for type 1 diabetes) for at least 15 minutes and ends when the glucose value is greater than or equal to the low limit of the target range for at least 15 minutes.
6. Serious Hypoglycaemia Event: An event starts when the glucose value remains below the hypoglycemia threshold (or TBR2 boundary value, e.g., 3.0 mmol/L or 54 mg/dL for type 1 diabetes) for at least 15 minutes, and ends when the glucose value stays at or above the hypoglycemia threshold (or TBR2 boundary value) for at least 15 minutes.
7. Hyperglycaemia Event: Event starts when the glucose value remains above the high limit of the target range (e.g., 10.0 mmol/L or 180 mg/dL for type 1 diabetes) for at least 15 minutes, and ends when the glucose value stays at or below the high limit of the target range for at least 15 minutes.
8. Serious Hyperglycaemia Event: Event starts when the glucose value remains above the high glucose threshold (or TAR2 boundary value, e.g., 13.9 mmol/L or 250 mg/dL for type 1 diabetes) for at least 15 minutes, and ends when the glucose value stays at or below the high glucose threshold (or TAR2 boundary value) for at least 15 minutes.
9. Highest-PG: Refers to peak glucose within 3 hours of a meal.
10. Tpeak: Refers to the time from mealtime to peak postprandial glucose.
11. PPGE: Equal to the postprandial glucose peak minus the preprandial glucose value.
12. Largest Amplitude of Glycemic Excursion(LAGE): Refers to the difference between the maximum and minimum glucose values during the monitoring period and evaluates the maximum glucose fluctuation.
13. Mean Amplitude of Glycemic Excursion(MAGE): After removing all glucose fluctuations whose amplitude did not exceed a certain threshold (1SDBG), the average value obtained by calculating the fluctuation amplitude based on the direction of the first valid fluctuation was used to evaluate the degree of intraday glucose fluctuation.
14. Mean of Daily Differences(MODD): It refers to the average of the absolute value of the difference between the relative measurements subtracted from each other over 2 consecutive days (i.e., a complete 48-hour period), reflecting the reproducibility of glucose from one day to another, and is used to evaluate the degree of fluctuation of glucose from day to day; MODD on the corresponding date reflects the degree of fluctuation of glucose values of the same day and the previous day.
15. Standard Deviation of Blood Glucose(SDBG): Refers to the standard deviation of measured glucose values during the monitoring date and evaluates the degree of overall deviation from the mean glucose value.
16. 50% median line: The median value of all glucose values over a given time period.
17. 25%-75% range(IQR): Represented by the 25th and 75th percentile curves, this area represents 50% of the glucose readings at any point in time and evaluates day-to-day glucose fluctuations, with wider intervals suggesting greater glucose fluctuations over the corresponding time period. The difference between the 25th and 75th percentile curves is shown here.
18. 5%-95% Range (IVR): Represented by the 5th and 95th percentile curves. This range represents 90% of glucose readings at any point in time, used to assess daily glucose fluctuation. A wider range indicates greater glucose fluctuation, often associated with irregular lifestyle patterns. This value represents the difference between the 5th and 95th percentile curves.
19. Ambulatory Glucose Profile (AGP): The AGP graph reflects your glucose fluctuations. The 50% median line shows the median of all glucose values at each time point — the smoother the curve, the more stable the glucose level. The 25%-75% shaded area covers 50% of glucose readings, and the 5%-95% shaded area covers 90% of glucose readings, indicating daily glucose fluctuations. A wider band represents greater fluctuations, while a narrower band indicates more stable glucose levels.

References

- [1] Clinical Targets for Continuous Glucose Monitoring Data Interpretation: Recommendations From the International Consensus on Time in Rang [J]. Diabetes Care 2019;42(8):1593–1603. <https://doi.org/10.2337/dci19-0028>
- [2] Glycemic Targets: Standards of Care in Diabetes—2023 [J].Diabetes Care 2023;46(Supplement_1):S97–S110.<https://doi.org/10.2337/dc23-S006>
- [3] International Diabetes Center captiurAGP® Version 5.0 [J] Final 3.11.2023
- [4] Continuous glucose monitoring and metrics for clinical trials: an international consensus statement [J].Lancet Diabetes Endocrinol 2023; 11: 42–57. [https://doi.org/10.1016/S2213-8587\(22\)00319-9](https://doi.org/10.1016/S2213-8587(22)00319-9)
- [5] Predictive models of post-prandial glucose response in persons with prediabetes and early onset type 2 diabetes: A pilot study [J]. Diabetes, Obesity and Metabolism 2025; 27(3):1515–1525. <https://doi.org/10.1111/dom.16160>

AGP Report: Continuous Glucose Monitoring




Name: *DEMO REPORT*

Reporting period: 15/09/2026 ~ 23/09/2026

Monitoring Days: 8 Days


Time CGM Active: 100%

Glucose Metrics	
Mean Glucose (MG)	113 mg/dL
Goal: <154mg/dL	
Glucose Management Indicator (GMI)	6 %
Goal: <7%	
Coefficient of variation (CV)	23.7 %
Defined as percent coefficient of variation	
Goal: ≤36%	



Daily Glucose Profiles

Each daily profile represents a midnight-to-midnight period.



ARCHPLUS HEALTHCARE PVT. LTD.

SHOP.47&47 A, G.H.B. BANDHAN COMPLEX,NR. ANKUR CHAR RASTA,NARANPURA, A'BAD

Patient: DEMO REPORT

Refd. By:

Pred.Eqns: RECORDERS

Date : 23-Dec-2023 09:56 PM

Age : 30 Years

Height : 162 Cms

Weight : 76 Kgs





ID: 20230818

Gender : Female

Smoker : Yes

Eth. Corr: 100

Temp : 98°C



Spirometry Results

Parameter	Pred	M.Pre	%Pred	M.Post	%Pred	%Imp
FVC (L)	02.60	03.22	124	03.32	128	+03
FEV1 (L)	02.08	02.56	123	02.76	133	+08
FEV1/FVC (%)	80.00	79.50	099	83.13	104	+05
FEF25-75 (L/s)	02.98	02.31	078	02.83	095	+23
PEFR (L/s)	06.63	06.40	097	06.89	104	+08
FIVC (L)	-----	02.86	---	02.77	---	-03
FEV.5 (L)	-----	01.95	---	02.13	---	+09
FEV3 (L)	02.53	03.09	122	03.22	127	+04
PIFR (L/s)	-----	02.26	---	02.41	---	+07
FEF75-85 (L/s)	-----	00.65	---	00.79	---	+22
FEF.2-1.2 (L/s)	05.40	05.12	095	06.14	114	+20
FEF 25% (L/s)	05.95	05.04	085	06.48	109	+29
FEF 50% (L/s)	04.66	02.68	058	03.19	068	+19
FEF 75% (L/s)	02.63	01.01	038	01.24	047	+23
FEV.5/FVC (%)	-----	60.56	---	64.16	---	+06
FEV3/FVC (%)	97.31	95.96	099	96.99	100	+01
FET (Sec)	-----	06.48	---	09.76	---	---
ExptTime (Sec)	-----	00.06	---	00.05	---	---
Lung Age (Yrs)	030	023	077	020	067	-13
FEV6 (L)	02.60	03.22	124	03.29	127	+02
FIF 25% (L/s)	-----	01.65	---	00.91	---	-45
FIF 50% (L/s)	-----	02.11	---	02.32	---	+10
FIF 75% (L/s)	-----	01.91	---	02.00	---	+05
SVC (L)	01.25	-----	---	-----	---	---
ERV (L)	-----	-----	---	-----	---	---
IRV (L)	-----	-----	---	-----	---	---
VE (L/min)	-----	-----	---	-----	---	---
Rf (l/min)	-----	-----	---	-----	---	---
Ti (sec)	-----	-----	---	-----	---	---
Te (sec)	-----	-----	---	-----	---	---
VT (L)	-----	-----	---	-----	---	---
VT/Ti	-----	-----	---	-----	---	---
Ti/Ttot	-----	-----	---	-----	---	---
IC (L)	-----	-----	---	-----	---	---
MVV (L/min)	111	-----	---	-----	---	---
MRf (l/min)	-----	-----	---	-----	---	---
MVT (L)	-----	-----	---	-----	---	---

Pre Medication Report Indicates

Spirometry within normal limits as (FEV1/FVC)%Pred >95 and FVC%Pred >80.

Post Medication Report Indicates

Spirometry within normal limits as (FEV1/FVC)%Pred >95 and FVC%Pred>80.

DR.